



## User's Manual

# Mouse IL-23R ELISA Kit

**REF** DEIABL646

 96T

**RUO**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of Mouse IL-23R.

### Reagents And Materials Provided

- 96-well plate precoated with anti-Mouse IL23r antibody: 1
- lyophilized recombinant Mouse IL23r standard: 10ng/tube
- biotinylated anti-Mouse IL23r antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

### Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

### Detection Range

78pg/ml-5000pg/ml

### Sensitivity

<10pg/ml