



User's Manual

Monkey GDF5/Bmp 14 ELISA Kit



DEIABL663



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of monkey primate GDF5.

Reagents And Materials Provided

- 96-well plate precoated with anti-Monkey Gdf5 antibody: 1
- lyophilized recombinant Monkey Gdf5 standard: 10ng/tube
- biotinylated anti-Monkey Gdf5 antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

Detection Range

31.2pg/ml-2000pg/ml

Sensitivity

<10pg/ml