



User's Manual

Human Plexin B3 ELISA Kit



DEIABL682



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of human Plexin B3.

Reagents And Materials Provided

- 96-well plate precoated with anti-Human PLXNB3 antibody: 1
- lyophilized recombinant Human PLXNB3 standard: 10ng/tube
- biotinylated anti-Human PLXNB3 antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

Detection Range

93.6pg/ml-6000pg/ml

Sensitivity

<15pg/ml



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com