



User's Manual

Mouse soluble CD6 ELISA Kit

REF

DEIABL706



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Sandwich High Sensitivity ELISA kit for Quantitative Detection of Mouse CD6.

Reagents And Materials Provided

- 96-well plate precoated with anti-Mouse Cd6 antibody: 1
- lyophilized recombinant Mouse Cd6 standard: 10ng/tube
- biotinylated anti-Mouse Cd6 antibody: 130ul (dilution 1:100)
- Avidin-Biotin-Peroxidase Complex(ABC): 130ul (dilution 1:100)
- sample diluent buffer: 30ml
- antibody diluent buffer: 12ml
- ABC diluent buffer: 12ml
- TMB color developing agent: 10ml
- TMB stop solution: 10ml
- Adhesive cover: 4

Storage

Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles(Shipped with wet ice.)

Detection Range

62.5pg/ml-4000pg/ml

Sensitivity

<10pg/ml



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com