



## User's Manual

# Thiamphenicol drug residue ELISA Kit



DEIANJ33



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Thiamphenicol drug residue ELISA Kit is intended for the quantitative determination of Thiamphenicol drug residue in raw milk, egg, seafood, tissue.

### Reagents And Materials Provided

1. Microtiter well plate: 8 microwells/strip, 12 strips/plate.
2. Standards
3. Washing Solution
4. Enzyme Conjugate Diluent
5. Substrate Solution (TMB)
6. Stop Solution

### Storage

Store all reagents at 2-8°C. Do not freeze. Reagents may be used until the expiration date on the box.