



User's Manual

Human Fluorescent Brightener KSN ELISA Kit

REF

DEIANJ34



96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Fluorescent Brightener KSN ELISA Kit is a competitive enzyme immunoassay for the quantitative analysis of Fluorescent Brightener KSN in serum, plasma and urine.

Reagents And Materials Provided

1. Microtiter plate (12 X 8 strips) coated with antibody, in an resealable aluminum pouch
2. Calibrators/Standards
3. HRP Conjugate
4. Antibody Solution
5. Wash Solution
6. Substrate (Color) Solution (TMB)
7. Stop Solution
8. Sample Diluent

Storage

Store the kit at 2 - 8°C.