



## User's Manual

# Butylated hydroxytoluene ELISA Kit



DEIANJ38



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Butylated hydroxytoluene ELISA Kit is a competitive enzyme immunoassay for the quantitative analysis of Butylated hydroxytoluene in oil.

### Reagents And Materials Provided

1. Microtiter plate (12 X 8 strips) coated with antibody, in an resealable aluminum pouch
2. Calibrators/Standards
3. HRP Conjugate
4. Antibody Solution
5. Wash Solution
6. Substrate (Color) Solution (TMB)
7. Stop Solution
8. Sample Diluent

### Storage

Store the kit at 2 - 8°C.

