



User's Manual

Hexadecadrol ELISA Kit



DEIANJ46



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Hexadecadrol ELISA Kit is a competitive enzyme immunoassay for the quantitative analysis of hexadecadrol in serum, plasma and urine.

Reagents And Materials Provided

1. Microtiter plate (12 X 8 strips) coated with antibody, in an resealable aluminum pouch
2. Calibrators/Standards
3. HRP Conjugate
4. Antibody Solution
5. Wash Solution
6. Substrate (Color) Solution (TMB)
7. Stop Solution
8. Sample Diluent

Storage

Store the kit at 2 - 8°C.

