



User's Manual

T-2 Toxin ELISA Kit



DEIAT2001



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This kit can be used in quantitative and qualitative analysis of T-2 toxin in feed.

Reagents And Materials Provided

1. Microtiter plate with 96 wells coated with antigen
2. Standard solutions(6 bottles:1ml/bottle): 0ppb, 1ppb, 3ppb, 9ppb, 27ppb, 81ppb
3. Concentrated enzyme conjugate 1.2ml
4. Enzyme conjugate diluent 12ml
5. Substrate solution A 7ml
6. Substrate Solution B 7ml
7. Stop solution 7ml

Storage

Storage condition: 2-8°C.

Storage period: 12months.

Precision

Variation coefficient of the ELISA kit is less than 10%.

Detection Limit

10ppb

Sensitivity

1ppb

Specificity

T-2 toxin: 100%