



User's Manual

Angiogenin ELISA Kit



DIA-XYA130



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION


Intended Use

The Human Angiogenin ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human Angiogenin concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Typical Standard Curve



Detection Range

16-500 pg/mL

Sensitivity

The Human Angiogenin ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human Angiogenin proteins within the range of 16-500 pg/mL.

Specificity

The Human Angiogenin ELISA is capable of recognizing both recombinant and naturally produced Human

Angiogenin proteins. Human Angiogenin was assayed against antigens at 50 ng/mL and did not exhibit significant cross-reactivity or interference.