



## User's Manual

# IGFBP-3 ELISA Kit



**DIA-XYA155**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION


### Intended Use

The Human IGFBP-3 ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human IGFBP-3 concentrations within any experimental sample including cell lysates, serum and plasma.

### Storage

4°C/6 Months

### Typical Standard Curve



### Detection Range

125-8000 pg/mL

### Sensitivity

The Human IGFBP-3 ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human IGFBP-3 proteins within the range of 125-8000 pg/mL.

### Specificity

The Human IGFBP-3 ELISA is capable of recognizing both recombinant and naturally produced Human IGFBP-3 proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross-reactivity or interference.

Human: IGF-I R, IGFBP-1, IGFBP-2, IGFBP-4, IGFBP-5, IGFBP-6