



User's Manual

MIP-3beta ELISA Kit

REF

DIA-XYA168



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION


Intended Use

The Human MIP-3 β ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human MIP-3 β concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Typical Standard Curve



Detection Range

16-1000 pg/mL

Sensitivity

The Human MIP-3 β ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human MIP-3 β proteins within the range of 16-1000 pg/mL.

Specificity

The Human MIP-3 β ELISA is capable of recognizing both recombinant and naturally produced Human MIP-3 β proteins. The antigens listed below were tested at 50 ng/mL and did not exhibit significant cross-reactivity or interference.

Human: MIP-3 β

Mouse: 6Ckine, MIP-1 α , MIP-1 β , TARC, TECK, RANTES