



User's Manual

CD163 ELISA Kit



DIA-XYA172



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION


Intended Use

The Human CD163 ELISA Kit contains the components necessary for quantitative determination of natural or recombinant Human CD163 concentrations within any experimental sample including cell lysates, serum and plasma.

Storage

4°C/6 Months

Typical Standard Curve



Detection Range

157-10000 pg/mL

Sensitivity

The Human CD163 ELISA Kit allows for the detection and quantification of endogenous levels of natural and/or recombinant Human CD163 proteins within the range of 157-10000 pg/mL.

Specificity

The Human CD163 ELISA is capable of recognizing both recombinant and naturally produced Human CD163

proteins. The antigens listed below were tested at 100 ng/mL and did not exhibit significant cross-reactivity or interference.

Human: CF6, TWEAK

Mouse: CF6, CD163, TWEAK