



Florfenicol Standard (99%) (DISNJ22)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	99%
Conjugate	N/A
Cas.No	73231-34-2
Applications	ELISA
Molecular Weight	358.21 g/mol
Format	Solid
Size	1 ml
Preservative	None
Storage	Room temperature

BACKGROUND

Introduction	Florfenicol is a fluorinated synthetic analog of thiamphenicol, mainly used in veterinary medicine.
Keywords	Florfenicol; Bovine respiratory disease; Aquaculture