



Chloramphenicol Standard (98%) (DISNJ29)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	98%
Conjugate	N/A
Cas.No	56-75-7
Applications	ELISA
Molecular Weight	323.13 g/mol
Format	Solid
Size	25g
Preservative	None
Storage	Room temperature

BACKGROUND

Introduction	Chloramphenicol is an antibiotic useful for the treatment of a number of bacterial infections.
Keywords	Chloremphenicol; chloramphenicol; Chloromycetin; Chlornitromycin; Levomycetin