



## Cefquinome Standard (98%) (DISNJ32)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Purity	98%
Conjugate	N/A
Cas.No	84957-30-2
Applications	ELISA
Molecular Weight	528.6 g/mol
Format	Solid
Size	5 mg
Preservative	None
Storage	Room temperature

### BACKGROUND

Introduction	Cefquinome is a fourth-generation cephalosporin with pharmacological and antibacterial properties valuable in the treatment of coliform mastitis and other infections.
Keywords	Cefquinome; Aminothiazolyl cephalosporin; PBPs