



Sulfamethazine Standard (98%) (DISNJ38)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	98%
Conjugate	N/A
Cas.No	57-68-1
Applications	ELISA
Molecular Weight	278.33 g/mol
Format	Solid
Size	25g
Preservative	None
Storage	Room temperature

BACKGROUND

Introduction	Sulfamethazine is a sulfonamide antibacterial.
Keywords	Sulfamethazine; Pyrimidine; Sulfanilamide