



Sulfamerazine Standard (98%) (DISNJ40)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	98%
Conjugate	N/A
Cas.No	127-79-7
Applications	ELISA
Molecular Weight	264.3 g/mol
Format	Solid
Size	250 mg
Preservative	None
Storage	Room temperature

BACKGROUND

Introduction	Sulfamerazine is a sulfonamide antibacterial.
Keywords	Sulfamerazine; Antibacterial agent; PABA