



# Thiamphenicol Standard (98%) (DISNJ45)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	98%
Conjugate	N/A
Cas.No	15318-45-3
Applications	ELISA
Molecular Weight	356.22 g/mol
Format	Solid
Size	1g
Preservative	None
Storage	Room temperature

## BACKGROUND

Introduction	Thiamphenicol is an antibiotic and immunosuppressive agent.
Keywords	Thiamphenicol; Monocarboxylic acid amide; Antimicrobial