



# Mouse anti-DENV NS1 Monoclonal Antibody (DMAB-A003)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Dengue NS1 antibody was raised in Mouse using native dengue virus NS1 as the immunogen.
<b>Specificity</b>	Specific for Dengue Virus Subtype 2
<b>Immunogen</b>	Native dengue virus NS1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Clone</b>	N64343
<b>Purification</b>	> 90% pure.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Lateral Flow
<b>Concentration</b>	Lot Specific
<b>Size</b>	500 µg
<b>Buffer</b>	Supplied in liquid form in 0.1M PBS buffer, pH 7.2, containing 0.09% sodium azide.
<b>Storage</b>	Store at 2-8 °C for short term storage. aliquot and store at -15 °C for long term storage, avoid repeated freeze/thaw cycles.

## BACKGROUND

**Keywords** Dengue Virus 1+2+3+4; Dengue Virus Type 1, 2, 3, 4; Group IV; Flaviviridae; Flavivirus;

Dengue virus; Dengue Virus; DENV;DENV NS1;Dengue Virus NS1;DENV Nonstructural Protein 1;Dengue Virus Nonstructural Protein 1;

---