



Mouse anti-Dengue NS1 Monoclonal Antibody (DMAB-A004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Dengue NS1 antibody was raised in Mouse using native dengue virus NS1 as the immunogen.
Specificity	Specific for Dengue Virus Subtype 2
Immunogen	Native dengue virus NS1
Isotype	IgG
Source/Host	Mouse
Clone	N64344
Purification	> 90% pure.
Conjugate	Unconjugated
Applications	Lateral Flow
Concentration	Lot Specific
Size	1 mg
Buffer	Supplied in liquid form in 0.1M PBS buffer, pH 7.2, containing 0.09% sodium azide.
Storage	Store at 2-8 °C for short term storage. aliquot and store at -15 °C for long term storage, avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Dengue Virus 1+2+3+4; Dengue Virus Type 1, 2, 3, 4; Group IV; Flaviviridae; Flavivirus;

Dengue virus; Dengue Virus; DENV;DENV NS1;Dengue Virus NS1;DENV Nonstructural Protein 1;Dengue Virus Nonstructural Protein 1;
