



Mouse anti-Ebola Virus Nucleoprotein monoclonal antibody (DMAB-A005)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Ebola Virus nucleoprotein.
Specificity	Reacts with the Ebola Virus nucleoprotein of Sudan and Zaire sequences. It reacts with Ebola Virus-Like Particles (VLPs).
Immunogen	Recombinant Ebola Virus nucleoprotein.
Source/Host	Mouse
Clone	C4263N
Purification	> 90% pure (SDS-PAGE). Protein A Chromatography.
Conjugate	Unconjugated
Applications	Suitable for use in ELISA, Lateral Flow and Western Blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot Specific
Size	1 mg
Buffer	0.01 M PBS, pH7.4
Preservative	None
Storage	Store at 2-8 °C.

BACKGROUND

Keywords

Ebola-Z nuclear protein; Ebola-Z; Nuclear protein; Ebola;Ebola Virus;Ebola Virus Nucleoprotein;Ebola Nucleoprotein;Ebola Virus NP
