



Mouse anti-Marburg Virus monoclonal antibody (DMAB-A008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG1
Source/Host	Mouse
Clone	IN835
Purification	>90%
Conjugate	Unconjugated
Applications	ELISA
Concentration	Lot Specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C.

BACKGROUND

Keywords	MBG; Marburg virus; Marburg; Marburgvirus; Mononegavirales; Filoviridae
----------	---