



Mouse anti-Yellow Fever NS1 monoclonal antibody (DMAB-A3009)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Culture Virus Infected Cell Lysate, Strain 17D.
Specificity	Yellow Fever Virus NS1 protein. No cross reaction with other Flavivirus proteins.
Immunogen	Culture Virus Infected Cell Lysate, Strain 17D
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	YFV
Clone	C3068N
Purification	> 90% pure (SDS-PAGE). Ion Exchange Chromatography.
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and Lateral Flow. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot Specific
Size	1 mg
Buffer	0.015 M Potassium Phosphate, 0.85% Sodium Chloride, pH 7.2.
Preservative	0.05% Sodium Azide

Storage

Store at 2-8 °C.

BACKGROUND

Keywords

Yellow Fever Virus NS1 antigen;YFV NS1 antigen;YFV NS1 protein;YFV Nonstructural protein;Yellow Fever Virus Nonstructural protein;Yellow Fever Virus;YFV
