



# Hi-Puri™ Rabbit Anti-Human ABCG8 Monoclonal antibody, clone 2E10 (DMAB- CDB25541)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	2E10 predominantly binds to the NBD of ABCG8. 2E10 inhibits ATP hydrolysis, likely by restricting the relative movement between the RecA and helical domain of ABCG8 NBD. mAb 2E10 inhibits the ATPase activity with an IC50 of 49.4 nM.
<b>Specificity</b>	Fab 2E10 interacts with both the RecA and the helical domains of the NBD from ABCG8.
<b>Target</b>	Human ABCG8
<b>Immunogen</b>	Full-length ABCG5/G8
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	2E10
<b>Purification</b>	>90% determined by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for use in SPR, ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid

<b>Concentration</b>	lot specific
<b>Size</b>	200 µg, 1 mg
<b>Buffer</b>	PBS (endotoxin < 1EU/mg, lower endotoxin levels may also be offered upon request)
<b>Preservative</b>	None
<b>Storage</b>	Short term at 2-8°C; long term storage in aliquots at -20°C; avoid freeze/thaw cycles.
<b>Ship</b>	Dry ice

## BACKGROUND

### Introduction

The ATP-binding cassette (ABC) transporters utilize the energy from ATP hydrolysis to transport a variety of substrates across membranes and are found in all kingdoms of life. In humans, there are 48 ABC transporters divided into 7 sub-families (A–G). While all ABC transporters consist of a pair of nucleotide-binding domains (NBDs) and a pair of transmembrane domains (TMDs), the subfamily G (ABCG) members possess a unique architecture in which the NBD is N-terminal to the TMD. Among the ABCG family members, ABCG5/G8 is essential for pumping cholesterol and phytosterols outward across apical membranes of enterocytes and hepatocytes. It serves as a gatekeeper of sterol transport and acts opposite to the Niemann–Pick C1-Like protein 1 (NPC1L1). NPC1L1 facilitates sterol influx across apical membranes of enterocytes from intestinal lumen and hepatocytes from bile canaliculus in the liver, whereas ABCG5/G8 mediates efflux of the sterols in the cells to allow proper absorption of cholesterol while restricting the absorption of structurally similar phytosterols. People with loss-of-function variants of ABCG5 or ABCG8 develop sitosterolemia, an autosomal disease characterized by impaired ability to eliminate dietary sterols. Sitosterolemic patients have considerably higher plasma levels of phytosterols, which can develop tendon xanthomas, and pose a high risk of cardiovascular disease. On the other hand, gain-of-function mutation variants of ABCG5/G8 are associated with gallstone disease.

### Keywords

ATP-binding cassette sub-family G member 8; ABCG8; STSL; sterolin 2; sterolin-2