



Hi-Puri™ Llama Anti-Human IL2RB Monoclonal antibody, clone Nb6 (DMAB-CDB26188)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Nb6 antibody specifically binds to the human IL-2R β ECD.
Target	Human IL2RB
Immunogen	The immunogen for Nb6 antibody is the human IL-2R β ECD (extracellular domain)
Isotype	VHH
Source/Host	Llama
Species Reactivity	Human
Clone	Nb6
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	Suitable for use in SPR, ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	lot specific
Size	200 μ g, 1 mg
Buffer	PBS (endotoxin < 1EU/mg, lower endotoxin levels may also be offered upon request)

Preservative	None
Storage	Short term at 2-8°C; long term storage in aliquots at -20°C; avoid freeze/thaw cycles.
Ship	Dry ice

BACKGROUND

Introduction IL2RB (Interleukin 2 Receptor Subunit Beta) is a Protein Coding gene. Diseases associated with IL2RB include Immunodeficiency 63 With Lymphoproliferation And Autoimmunity and Ichthyosis. Among its related pathways are TGF-Beta Pathway and MIF Mediated Glucocorticoid Regulation.

Keywords CD122; IL2RB; Interleukin 2 Receptor Subunit Beta