



Mouse Anti-HPV L1 Monoclonal antibody, clone 56B5 (DMAB-CS23002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	HPV L1
Target	HPV L1
Immunogen	Recombinant HPV L1 VLP from E. coli
Isotype	IgG
Source/Host	Mouse
Species Reactivity	HPV
Clone	56B5
Purification	Protein G
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	Tris-Gly
Preservative	None
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction

Human papillomaviruses from the family Papillomaviridae belong to a category of small non-enveloped circular double-stranded DNA viruses, measuring 50–55 nm in diameter. The late gene coding section has two parts: L1 and L2. L1 codes for a major viral capsid protein, while L2 codes for a minor viral capsid structure.

Keywords

Human Papillomavirus; HPV; HPV L1
