



# Magic™ Mouse Anti-Human TFF2 Monoclonal antibody, clone TTF1 (DMAB-CS24605)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Human TFF2
<b>Target</b>	Human TFF2
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	TTF1
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(Cap), LFIA, CLIA, TINIA We recommend the following antibodies for sandwich immunoassay (Capture - Detection): DMAB-CS24605 - DMAB-CS24606
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	Supplied as a 0.2 µM filtered solution of PBS, PH7.4.
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

## BACKGROUND

<b>Introduction</b>	TFF2 (Trefoil Factor 2) is a Protein Coding gene. Diseases associated with TFF2 include Gastritis and Atrophic Gastritis. Among its related pathways are Pancreatic cancer subtypes and Gastrin signaling pathway. Gene Ontology (GO) annotations related to this gene include CXCR4 chemokine receptor binding. An important paralog of this gene is TFF1.
<b>Keywords</b>	TFF2; trefoil factor 2; SP; SML1; spasmodysin; spasmodytic protein 1; spasmodytic polypeptide