



Magic™ Mouse Anti-Human AMBP Monoclonal antibody, clone B2N1 (DMAB-CS24613)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human AMBP
Target	Human AMBP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	B2N1
Conjugate	Unconjugated
Applications	ELISA(Cap), LFIA, CLIA, TINIA, POCT We recommend the following antibodies for sandwich immunoassay (Capture - Detection): DMAB-CS24613 - DMAB-CS24614
Format	Liquid
Size	1 mg
Buffer	Supplied as a 0.2 µM filtered solution of PBS, PH7.4.
Preservative	None
Storage	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

BACKGROUND

Introduction

AMBP (Alpha-1-Microglobulin/Bikunin Precursor) is a Protein Coding gene. Diseases associated with AMBP include Hepatitis E and Obstructive Nephropathy. Among its related pathways are Binding and Uptake of Ligands by Scavenger Receptors and Vesicle-mediated transport. Gene Ontology (GO) annotations related to this gene include protein homodimerization activity and serine-type endopeptidase inhibitor activity. An important paralog of this gene is WFIKK2.

Keywords

Alpha-1-microglobulin; AMBP; A1M; HCP; ITI; UTI; EDC1; HI30; ITIL; IATIL; ITILC
