



# Magic™ Human Anti-Human VSNL1

## Monoclonal antibody, clone WMQ3 (DMAB-CS24625)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	Human VSNL1
<b>Target</b>	Human VSNL1
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	WMQ3
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(Cap), LFIA, CLIA, TINIA, POCT We recommend the following antibodies for sandwich immunoassay (Capture - Detection): DMAB-CS24625 - DMAB-CS24626
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	Supplied as a 0.2 µM filtered solution of PBS, PH7.4.
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

# BACKGROUND

Introduction	VSNL1 (Visinin Like 1) is a Protein Coding gene. Diseases associated with VSNL1 include Abcd Syndrome. Among its related pathways is Neuroscience. Gene Ontology (GO) annotations related to this gene include calcium ion binding. An important paralog of this gene is HPCAL4.
Keywords	VSNL1; visinin-like 1; HLP3; VILIP; HPCAL3; HUVISL1; VILIP-1; visinin-like protein 1; VLP-1; hippocalcin-like protein 3