



HiResNb™ Alpaca Anti-Human CCL5 Monoclonal Antibody, clone A5043 (DMAB- CS24770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human CCL5/RANTES
Target	Human CCL5
Isotype	VHH
Source/Host	Alpaca
Species Reactivity	Human
Clone	A5043
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA, IP
Format	Liquid
Size	100 µg, 1 mg
Buffer	10mM PBS (pH 7.5)
Preservative	None
Storage	Store at -20°C (Avoid freeze / thaw cycles).

BACKGROUND

Introduction

Chemokine (C-C motif) ligand 5 (also CCL5) is a protein which in humans is encoded by the CCL5 gene. It is also known as RANTES (regulated on activation, normal T cell expressed and secreted).

Keywords

CCL5; chemokine (C-C motif) ligand 5; SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta; C-C motif chemokine 5
