



Magic™ Mouse Anti-Canine Progesterone Monoclonal antibody, clone DQ02 (DMAB- CS25011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Canine Progesterone
Target	Canine Progesterone
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Canine
Clone	DQ02
Conjugate	Unconjugated
Applications	ELISA(Det), LFIA, CLIA, TINIA, POCT We recommend the following antibodies for sandwich immunoassay (Capture - Detection): DMAB-CS25010 - DMAB-CS25011
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from sterile PBS, PH 7.4
Preservative	None
Storage	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

BACKGROUND

Introduction

Progesterone also known as P4 (pregn-4-ene-3,20-dione) is a C-21 steroid hormone involved in the female menstrual cycle, pregnancy (supports gestation) and embryogenesis of humans and other species. Progesterone belongs to a class of hormones called progestogens, and is the major naturally occurring human progestogen.

Keywords

Progesterone; P4
