



Magic™ Mouse Anti-Feline CRP Monoclonal antibody, clone GD01 (DMAB-CS25018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Feline CRP
Target	Feline CRP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Feline
Clone	GD01
Conjugate	Unconjugated
Applications	ELISA(Cap), LFIA, CLIA, TINIA, POCT We recommend the following antibodies for sandwich immunoassay (Capture - Detection): DMAB-CS25018 - DMAB-CS25019
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from sterile PBS, PH 7.4
Preservative	None
Storage	Store at -20°C to -80°C under sterile conditions. Avoid repeated freeze-thaw cycles.

BACKGROUND

Introduction

Cat/Feline CRP (C-reactive protein) is a specific biomarker found in cats and feline species. It is produced by the liver and gets released into the bloodstream in response to inflammation. This biomarker plays a crucial role in monitoring and assessing inflammatory processes that occur within the bodies of these animals.

Keywords

CRP; C-reactive protein
