



Mouse Anti-EV71 monoclonal antibody, clone N2D21 (DMAB-EV71)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Isotype | IgG1 |
| Source/Host | Mouse |
| Species Reactivity | EV |
| Clone | N2D21 |
| Purification | > 90% pure (SDS-PAGE) |
| Conjugate | Unconjugated |
| Applications | ELISA, CLIA, LFIA |
| Molecular Weight | 150 KD |
| Format | Solution |
| Concentration | Lot Specific |
| Size | 1 mg |
| Buffer | 20 mM PBS (PH 7.4) |
| Preservative | None |
| Storage | Store at 2-8°C and at -20°C for long-term storage. Don't repeated freeze-thaw. |

BACKGROUND

Keywords

Enterovirus 71 MRS 3663149/2010;Picornaviridae;Enterovirus 71;EV71;EV 71;EV-71
