



Rabbit Anti-Human CCL18 Monoclonal Antibody, Clone S16-7C5 (DMAB-JX2361713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S16-7C5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Capture)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361713-DMAB-JX2361714</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Chemotactic factor that attracts lymphocytes but not monocytes or granulocytes. May be involved in B-cell migration into B-cell follicles in lymph nodes. Attracts naive T-lymphocytes toward dendritic cells and activated macrophages in lymph nodes, has chemotactic activity for naive T-cells, CD4+ and CD8+ T-cells and thus may play a role in both humoral and cell-mediated immunity responses.
---------------------	--

GENE INFORMATION

Entrez Gene ID	6362
-----------------------	----------------------

UniProt ID	P55774
-------------------	------------------------