



Rabbit Anti-Human KRT19 Monoclonal Antibody, Clone S18-6O2 (DMAB-JX2361717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S18-6O2
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA (Capture) We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361717-DMAB-JX2361718 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.
---------------------	--

GENE INFORMATION

Entrez Gene ID	3880
-----------------------	----------------------

UniProt ID	P08727
-------------------	------------------------