



Rat Anti-Human FGA Monoclonal Antibody, Clone 10B3 (DMAB-JX2361722)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthesized peptide derived from FGA
Isotype	IgG1
Source/Host	Rat
Species Reactivity	Human
Clone	10B3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Detection)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361721-DMAB-JX2361722</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Fibrinogen has a double function: yielding monomers that polymerize into fibrin and acting as a cofactor in platelet aggregation.
---------------------	---

GENE INFORMATION

Entrez Gene ID	2243
-----------------------	----------------------

UniProt ID	P02671
-------------------	------------------------