



Rabbit Anti-Human CXCL1 Monoclonal Antibody, Clone S19-9E1 (DMAB-JX2361728)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S19-9E1
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Detection)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361727-DMAB-JX2361728</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Has chemotactic activity for neutrophils. May play a role in inflammation and exerts its effects on endothelial cells in an autocrine fashion. In vitro, the processed forms GRO-alpha(4-73), GRO-alpha(5-73) and GRO-alpha(6-73) show a 30-fold higher chemotactic activity.
---------------------	---

GENE INFORMATION

Entrez Gene ID	2919
-----------------------	----------------------

UniProt ID	P09341
-------------------	------------------------