



Rabbit Anti-Mouse IL2 Monoclonal Antibody, Clone S17-9Y5 (DMAB-JX2361747)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S17-9Y5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA (Capture) We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361747-DMAB-JX2361748 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	IL-2 is a potent lymphoid cell growth factor which exerts its biological activity primarily on T cells. Additionally, IL-2 has been found to stimulate growth and differentiation of B cells, NK cells, LAK cells, monocytes, and oligodendrocytes.
---------------------	---

GENE INFORMATION

Entrez Gene ID	16183
-----------------------	-----------------------

UniProt ID	P04351
-------------------	------------------------