



Rabbit Anti-Mouse IL2 Monoclonal Antibody, Clone S14-7G6 (DMAB-JX2361748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S14-7G6
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Detection)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361747-DMAB-JX2361748</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	IL-2 is a potent lymphoid cell growth factor which exerts its biological activity primarily on T cells. Additionally, IL-2 has been found to stimulate growth and differentiation of B cells, NK cells, LAK cells, monocytes, and oligodendrocytes.
---------------------	---

GENE INFORMATION

Entrez Gene ID	16183
-----------------------	-----------------------

UniProt ID	P04351
-------------------	------------------------