



Rabbit Anti-Mouse IL8 Monoclonal Antibody, Clone S18-7B3 (DMAB-JX2361760)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S18-7B3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA (Detection) We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361759-DMAB-JX2361760 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an inflammatory stimulus. IL-8(6-77) has a 5-10-fold higher activity on neutrophil activation, IL-8(5-77) has increased activity on neutrophil activation and IL-8(7-77) has a higher affinity to receptors CXCR1 and CXCR2 as compared to IL-8(1-77), respectively.
---------------------	---

GENE INFORMATION

Entrez Gene ID	3576
-----------------------	----------------------

UniProt ID	P10145
-------------------	------------------------