



Rabbit Anti-Human MIF Monoclonal Antibody, Clone S18-2Y3 (DMAB-JX2361770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S18-2Y3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Detection)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361769 - DMAB-JX2361770</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	MIF promotes the systemic inflammatory response by counter-regulating glucocorticoid-mediated inhibition of immune-cell activation and proinflammatory cytokine production. MIF may mediate tissue destruction through the induction of proteinases.
---------------------	--

GENE INFORMATION

Entrez Gene ID	4282
-----------------------	----------------------

UniProt ID	P14174
-------------------	------------------------