



Rabbit Anti-Human CCL3 Monoclonal Antibody, Clone S15-0A5 (DMAB-JX2361772)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S15-0A5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Detection)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361771-DMAB-JX2361772</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Monokine with inflammatory and chemokinetic properties. Binds to CCR1, CCR4 and CCR5. One of the major HIV-suppressive factors produced by CD8+ T-cells. Recombinant MIP-1-alpha induces a dose-dependent inhibition of different strains of HIV-1, HIV-2, and simian immunodeficiency virus (SIV).
---------------------	---

GENE INFORMATION

Entrez Gene ID	6348
-----------------------	----------------------

UniProt ID	P10147
-------------------	------------------------