



# Rabbit Anti-Human CCL4 Monoclonal Antibody, Clone S17-3P7 (DMAB-JX2361773)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	S17-3P7
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	<p>ELISA (Capture)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361773-DMAB-JX2361774</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µg
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	None

<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

<b>Ship</b>	Wet ice
-------------	---------

---

## BACKGROUND

<b>Introduction</b>	Monokine with inflammatory and chemokinetic properties. Binds to CCR5. One of the major HIV-suppressive factors produced by CD8+ T-cells. Recombinant MIP-1-beta induces a dose-dependent inhibition of different strains of HIV-1, HIV-2, and simian immunodeficiency virus (SIV).
---------------------	---

---

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">6351</a>
-----------------------	----------------------

<b>UniProt ID</b>	<a href="#">P13236</a>
-------------------	------------------------