



Rabbit Anti-Human CCL5 Monoclonal Antibody, Clone S17-9W0 (DMAB-JX2361775)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S17-9W0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	<p>ELISA (Capture)</p> <p>We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361775-DMAB-JX2361776</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Chemoattractant for blood monocytes, memory T-helper cells and eosinophils. Causes the release of histamine from basophils and activates eosinophils. Binds to CCR1, CCR3, CCR4 and CCR5.
---------------------	---

GENE INFORMATION

Entrez Gene ID	6352
-----------------------	----------------------

UniProt ID	P13501
-------------------	------------------------