



Rabbit Anti-Mouse VEGFA Monoclonal Antibody, Clone S12-2M4 (DMAB-JX2361786)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	S12-2M4
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA (Detection) We recommend the following as antibody pair (Capture - Detection): DMAB-JX2361785-DMAB-JX2361786 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	0.01M PBS, pH 7.2
Preservative	None

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Growth factor active in angiogenesis, vasculogenesis and endothelial cell growth. Induces endothelial cell proliferation, promotes cell migration, inhibits apoptosis and induces permeabilization of blood vessels. Binds to the FLT1/VEGFR1 and KDR/VEGFR2 receptors, heparan sulfate and heparin. May play a role in increasing vascular permeability during lactation, when increased transport of molecules from the blood is required for efficient milk protein synthesis (By similarity). Binding to NRP1 receptor initiates a signaling pathway needed for motor neuron axon guidance and cell body migration, including for the caudal migration of facial motor neurons from rhombomere 4 to rhombomere 6 during embryonic development.
---------------------	--

GENE INFORMATION

Entrez Gene ID	22339
-----------------------	-----------------------

UniProt ID	Q00731
-------------------	------------------------