



Rabbit Anti-Human NCOA1 Monoclonal Antibody, Clone S18-6C8 (DMAB-JX2362001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	KAT13A
Immunogen	Recombinant protein of human SRC1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S18-6C8
Purification	Affinity Purification
Conjugate	Unconjugated
Applications	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	20 µL, 50 µL, 100 µL
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Preservative	0.01% Sodium azide

Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Introduction	Nuclear receptor coactivator that directly binds nuclear receptors and stimulates the transcriptional activities in a hormone-dependent fashion. Involved in the coactivation of different nuclear receptors, such as for steroids (PGR, GR and ER), retinoids (RXRs), thyroid hormone (TRs) and prostanoids (PPARs).
---------------------	---

GENE INFORMATION

Entrez Gene ID	8648
-----------------------	----------------------

UniProt ID	Q15788
-------------------	------------------------