



Magic™ Mouse anti-Human cardiac troponin I (cTnI) Monoclonal antibody, clone 2H22 (DMAB-L21061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with phosphorylated form of cTnI.
Antigen Description	Troponin I (TnI), along with troponin T (TnT) and troponin C (TnC), is one of 3 subunits that form the troponin complex of the thin filaments of striated muscle. TnI is the inhibitory subunit; blocking actin-myosin interactions and thereby mediating striated muscle relaxation. The TnI subfamily contains three genes: TnI-skeletal-fast-twitch, TnI-skeletal-slow-twitch, and TnI-cardiac. This gene encodes the TnI-cardiac protein and is exclusively expressed in cardiac muscle tissues. Mutations in this gene cause familial hypertrophic cardiomyopathy type 7 (CMH7) and familial restrictive cardiomyopathy (RCM). [provided by RefSeq, Jul 2008]
Specificity	Phosphorylated form of cTnI.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2H22
Purification	Chromatography on protein A Sepharose
Conjugate	Unconjugated
Applications	MAb is working in Western blotting and ELISA.
Buffer	PBS, pH 7.4, 0.1 % sodium azide (NaN ₃)

Preservative	0.1% Sodium Azide
--------------	-------------------

Storage	4°C
---------	-----

GENE INFORMATION

Gene Name	TNNI3 troponin I type 3(cardiac) [Homo sapiens]
-----------	---

Official Symbol	TNNI3
-----------------	-------

Synonyms	TNNI3; troponin I type 3 (cardiac); CMH7; RCM1; cTnI; CMD2A; TNNC1; MGC116817; troponin I, cardiac; Troponin I, cardiac muscle; Cardiac troponin I; familial hypertrophic cardiomyopathy 7
----------	--

Entrez Gene ID	7137
----------------	----------------------

Protein Refseq	NP_000354
----------------	---------------------------

UniProt ID	P19429
------------	------------------------

Chromosome Location	19q13.4
---------------------	---------

Pathway	Hypertrophic cardiomyopathy (HCM); Muscle contraction; Cardiac muscle contraction
---------	---

Function	actin binding; calcium-dependent protein binding; protein domain specific binding; protein kinase binding; troponin C binding; troponin T binding; calcium channel inhibitor activity
----------	---