



# Magic™ Anti-HPV type 16 Monoclonal antibody, Clone D6 (DMAB-L21074)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	HPV type 16
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	HPV
Clone	D6
Purification	Chromatography on protein G Sepharose
Conjugate	Unconjugated
Applications	Indirect EIA with recombinant protein.
Buffer	PBS, pH 7.4, 0.1 % sodium azide (NaN <sub>3</sub> )
Preservative	0.1% Sodium Azide
Storage	4°C

## BACKGROUND

Introduction	Human papillomavirus (HPV) is a virus from the papillomavirus family that is capable of infecting humans. Like all papillomaviruses, HPVs establish productive infections only in keratinocytes of the skin or mucous membranes. While the majority of the know
Keywords	Human papillomavirus; HPV;