



# Magic™ Anti-IFN- $\gamma$ Monoclonal antibody, Clone F4G11 (DMAB-S002)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Reacts with: African and Asian elephant IFN- $\gamma$ , cross-reacts with horse
Immunogen	recombinant elephant rEplFN- $\gamma$ (expressed in HEK293 cells)
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	elephant
Purification	GammaBindPlus purified
Conjugate	Unconjugated
Applications	ELISA(Dec), WB We recommend the following for sandwich ELISA (Capture - Detection): DMAB-S001 - DMAB-S002
Format	Liquid
Concentration	0.5 mg/ml
Size	100 $\mu$ g
Buffer	PBS
Preservative	None
Storage	Store at 4° C short term (1-2 weeks) Aliquot and store at -21° C

## GENE INFORMATION

Entrez Gene ID [100663285](#)

---

UniProt ID [G3SUC6](#)

---